



A laryngoscope is used to examine and visualise a patient's upper airway and aid placement of a tracheal tube during intubations. This product must only be used by a qualified physician. (N.B. USA and Canadian Federal Law restricts this device to sale by or on the order of a physician.)

- Callisto blades and handles conform to the ISO 7376 standards (fibre optic fitting).
- Use only with handles or blades that conform to ISO 7376 standards (fibre optic fitting).
- Supplied non-sterile, packed in clean room conditions.
- Do not reuse, sterilise or reprocess.

Pre-use Checks – before use on a patient

- Visually inspect the product for any sharp edges, burrs or other damage.
- After inserting batteries (see images for instructions) depress the bulb collar on the handle to confirm that the bulb illuminates
- Fit the blade to the handle and check the light functions correctly.

IMPORTANT NOTE: Always keep spare blades and handles available in case of failure or emergency.

Cautions

- This device is designed to be SINGLE-USE and should only be used on a single patient during a single procedure and then discarded. Do not reuse or reprocess.
- Bulbs, if left illuminated in an exposed position, may generate heat sufficient to burn human tissue.
- Do not use in the vicinity of intense magnetic fields such as MRI scanners.
- Power outputs from some rechargeable cells may fall rapidly during use, resulting in rapid failure of illumination.
- This product is designed to be disposable with a limited operational life. Repeated testing of this equipment before use could significantly reduce the life of the product or result in operational failure.

Fitting & Removing a Laryngoscope Blade to a Laryngoscope Handle

1. Hook the blade onto the hook bar on the handle (1)
2. Apply pressure to the heel of the blade until the blade 'clicks' into position over the bar (2)
3. Lift the tip of the blade until it locks into the operating position and the bulb illuminates (3)
4. To remove the blade, reverse the instructions above.

Cautionary Note: It is possible to incorrectly mount the blade. This may result in a non-functioning or damaged blade and/or handle.

Fitting Batteries

Use alkaline dry cell batteries. Use either 2 'C' size batteries for medium handles or 2 'AA' size for small handles. (The use of rechargeable batteries are not recommended in these handles.)

1. Unscrew the head block from the battery compartment (4).
2. Insert new batteries as shown (5).
3. Replace and securely tighten block.

Please ensure that used blades, handles and batteries are discarded as clinical waste in accordance with local policies and all relevant regulatory guidelines.

CODE	DESCRIPTION
DS.3940.150.05	Callisto Macintosh Blade No.0 (Box of 10)
DS.3940.150.10	Callisto Macintosh Blade No.1 (Box of 10)
DS.3940.150.15	Callisto Macintosh Blade No.2 (Box of 10)
DS.3940.150.20	Callisto Macintosh Blade No.3 (Box of 10)
DS.3940.150.25	Callisto Macintosh Blade No.4 (Box of 10)
DS.3940.185.05	Callisto Miller Blade No.0 (Box of 10)
DS.3940.185.10	Callisto Miller Blade No.1 (Box of 10)
DS.3940.185.15	Callisto Miller Blade No.2 (Box of 10)
DS.3940.185.20	Callisto Miller Blade No.3 (Box of 10)
DS.3940.235.10	Callisto Robertshaw Blade No.1 (Box of 10)
DS.3940.240.10	Callisto Seward Blade No.1 (Box of 10)
DS.3940.350.10	Callisto S Handle Medium 'C' (Box of 10)



Non-Sterile



Do not re-use



Do not re-sterilise



Do not use if packaging has been opened or damaged.



Latex Free



Federal (U.S.A) law restricts this device to sale by or on the order of a physician.



Consult instructions for use





CODE	DESCRIPTION
DS.3940.370.10	Callisto LED Handle Medium 'C' (Box of 10)

	Non-Sterile
	Do not re-use
	Do not re-sterilise
	Do not use if packaging has been opened or damaged.
	Latex Free
Rx ONLY	Federal (U.S.A) law restricts this device to sale by or on the order of a physician.
	Consult instructions for use
CE	

A laryngoscope is used to examine and visualise a patient's upper airway and aid placement of a tracheal tube during intubations. This product must only be used by a qualified physician. (N.B. USA and Canadian Federal Law restricts this device to sale by or on the order of a physician.)

- This Callisto LED handle conforms to the ISO 7376 standard (fibre optic fitting).
- Use only with blades that conform to ISO 7376 standard (fibre optic fitting).
- Supplied non-sterile, packed in clean room conditions.
- Do not reuse, sterilise or reprocess.

Pre-use Checks – before use on a patient

- Visually inspect the product for any sharp edges, burrs or other damage.
- After inserting batteries (see images for instructions) depress the bulb collar on the handle to confirm that the bulb illuminates
- Fit the blade to the handle and check the light functions correctly.

IMPORTANT NOTE: Always keep spare blades and handles available in case of failure or emergency.

Cautions

- This device is designed to be SINGLE-USE and should only be used on a single patient during a single procedure and then discarded. Do not reuse or reprocess.
- Do not use in the vicinity of intense magnetic fields such as MRI scanners.
- Power outputs from some rechargeable cells may fall rapidly during use, resulting in rapid failure of illumination.
- This product is designed to be disposable with a limited operational life. Repeated testing of this equipment before use could significantly reduce the life of the product or result in operational failure.
- Users should not look directly into the LED and avoid pointing the LED handle towards other staff or the patient.

Fitting & Removing a Laryngoscope Blade to a Laryngoscope Handle

1. Hook the blade onto the hook bar on the handle (1)
2. Apply pressure to the heel of the blade until the blade 'clicks' into position over the bar (2)
3. Lift the tip of the blade until it locks into the operating position and the bulb illuminates (3)
4. To remove the blade, reverse the instructions above

Cautionary Note: It is possible to incorrectly mount the blade. This may result in a non-functioning or damaged blade and or handle.

Inserting the batteries

Use 2 'C' size alkaline dry cell batteries. (The use of rechargeable batteries are not recommended in these handles)

1. Unscrew the end cap from the bottom of the handle (4)
2. Insert new batteries as shown (5)
3. Replace and securely tighten the end cap.

Please ensure that used blades, handles and batteries are discarded as clinical waste in accordance with local policies and all relevant regulatory guidelines.